

# Clinical Data Visualizer – User Guide

Clinical Data Visualizer Team

February 18, 2026

## Contents

|           |                                                         |           |
|-----------|---------------------------------------------------------|-----------|
| <b>1</b>  | <b>Introduction</b>                                     | <b>3</b>  |
| 1.1       | Key Features . . . . .                                  | 3         |
| 1.2       | Supported Data Sources . . . . .                        | 3         |
| <b>2</b>  | <b>Launching the Application</b>                        | <b>4</b>  |
| 2.1       | Starting the App . . . . .                              | 4         |
| 2.2       | Application Overview . . . . .                          | 4         |
| <b>3</b>  | <b>Preparing Your Data</b>                              | <b>5</b>  |
| 3.1       | Patient Folder Structure . . . . .                      | 5         |
| 3.2       | Folder Naming Rules . . . . .                           | 5         |
| 3.3       | Expected File Types . . . . .                           | 5         |
| <b>4</b>  | <b>Loading Database Options</b>                         | <b>7</b>  |
| 4.1       | Option 1: Default Visualization (Quick Start) . . . . . | 7         |
| 4.2       | Option 2: Custom Configuration File . . . . .           | 7         |
| <b>5</b>  | <b>Configuring Patient Options</b>                      | <b>8</b>  |
| 5.1       | Global Options . . . . .                                | 8         |
| 5.2       | Per-Source Options . . . . .                            | 8         |
| <b>6</b>  | <b>Processing the Visualization</b>                     | <b>9</b>  |
| 6.1       | What Happens During Processing . . . . .                | 9         |
| <b>7</b>  | <b>Interacting with Plots</b>                           | <b>10</b> |
| 7.1       | Navigation Controls . . . . .                           | 10        |
| 7.2       | Zooming and Panning . . . . .                           | 10        |
| 7.3       | Dynamic Resampling (FigureResampler) . . . . .          | 10        |
| <b>8</b>  | <b>Annotations and Shape Tools</b>                      | <b>11</b> |
| 8.1       | Drawing Annotations . . . . .                           | 11        |
| 8.2       | Shape Management . . . . .                              | 11        |
| 8.3       | Annotation Properties . . . . .                         | 11        |
| 8.4       | Persistence . . . . .                                   | 11        |
| <b>9</b>  | <b>Configuration File Reference</b>                     | <b>12</b> |
| 9.1       | patient_options.json . . . . .                          | 12        |
| 9.2       | database_options.json . . . . .                         | 12        |
| <b>10</b> | <b>Troubleshooting</b>                                  | <b>14</b> |
| 10.1      | Browser Does Not Open Automatically . . . . .           | 14        |
| 10.2      | No Data Found . . . . .                                 | 14        |

|                                              |           |
|----------------------------------------------|-----------|
| 10.3 Slow Loading . . . . .                  | 14        |
| 10.4 Time Alignment Issues . . . . .         | 14        |
| 10.5 Application Crashes or Errors . . . . . | 14        |
| <b>11 Appendix: Supported Data Sources</b>   | <b>15</b> |

# 1 Introduction

Clinical Data Visualizer is an interactive dashboard for visualizing clinical physiological signals. It allows clinicians and researchers to explore, compare, and annotate time-series data from multiple medical devices in a single unified interface.

## 1.1 Key Features

- **Multi-source visualization:** Display signals from up to 9 different clinical data sources simultaneously.
- **Interactive plots:** Zoom, pan, and explore data at any time scale with automatic resampling.
- **Annotations:** Draw lines and rectangles directly on plots to mark events or regions of interest. Annotations are saved and persist across sessions.
- **Flexible configuration:** Choose which signals to display, customize labels, units, colors, and group related signals together.
- **Export:** Generate standalone HTML visualizations for sharing.

## 1.2 Supported Data Sources

| Data Source        | Description                                                           |
|--------------------|-----------------------------------------------------------------------|
| Philips Waves      | High-frequency waveform data from Philips patient monitors            |
| Philips Numerics   | Numeric parameter data from Philips monitors (heart rate, SpO2, etc.) |
| EIT (PulmoVista)   | Electrical Impedance Tomography data from Draeger PulmoVista          |
| FluxMed Signals    | Waveform data from FluxMed respiratory monitors                       |
| FluxMed Parameters | Parameter data from FluxMed monitors                                  |
| Servo-U            | Ventilator data from Getinge Servo-U                                  |
| Mindray            | Patient monitor data from Mindray scopes                              |
| Syringe            | Syringe pump infusion data                                            |
| Other (Generic)    | Auto-discovers CSV or Parquet files with datetime columns             |

## 2 Launching the Application

### 2.1 Starting the App

Locate the **ClinicalVisuAppAlexis** executable in the application folder and double-click it.

A terminal window will appear showing the application starting up. After a few seconds, your default web browser will automatically open at:

`http://127.0.0.1:8050`

If the browser does not open automatically, manually navigate to the address above.

Application launch screen

Figure 1: Application launch screen

### 2.2 Application Overview

The interface is organized top-to-bottom in the following order:

1. **Database Options** – Load or select a visualization configuration.
2. **Patient Options** – Configure data folder, time range, and per-source settings.
3. **Process Button** – Start the visualization.
4. **Shape Controls** – Manage annotations (visible after processing).
5. **Visualization Area** – Interactive plots.

Application main interface

Figure 2: Application main interface

## 3 Preparing Your Data

### 3.1 Patient Folder Structure

Each patient’s data must be organized in a root folder with **one subfolder per data source**. The application automatically identifies data sources based on keywords in subfolder names.

```
Patient1/
  philips_waves/      Philips waveform data (.parquet)
  philips_numerics/   Philips numeric data
  eit/                EIT PulmoVista data (.asc)
  fluxmed_signals/    FluxMed waveform data
  fluxmed_parameters/ FluxMed parameter data
  servo_u/            Servo-U ventilator data (.sta)
  mindray/            Mindray scope data (.xml or .csv)
  syringe/            Syringe pump data
  other/              Generic data (.csv or .parquet)
  tdv_visu/           Auto-created: cached data and outputs
```

You only need subfolders for the data sources you actually have. Empty or missing subfolders are simply skipped.

### 3.2 Folder Naming Rules

Folder names are **flexible** – they just need to contain the required keywords (case-insensitive, any separator allowed):

| Data Source        | Required Keywords    | Examples                                    |
|--------------------|----------------------|---------------------------------------------|
| Philips Waves      | philips + waves      | philips_waves, Philips-Waves, waves_philips |
| Philips Numerics   | philips + numerics   | philips_numerics, Philips-Numerics          |
| EIT                | eit                  | eit, EIT, EIT_Data                          |
| FluxMed Signals    | fluxmed + signals    | fluxmed_signals, FluxMed-Signals            |
| FluxMed Parameters | fluxmed + parameters | fluxmed_parameters, FluxMed_Parameters      |
| Servo-U            | servo                | servo_u, Servo-U, SERV0                     |
| Mindray            | mindray              | mindray, Mindray                            |
| Syringe            | syringe              | syringe, Syringe, syringe_pumps             |
| Other              | other                | other, Other                                |

### 3.3 Expected File Types

- **Philips Waves:** .parquet files
- **Philips Numerics:** Files containing “numerics” in the filename
- **EIT:** .asc files
- **FluxMed:** Files containing “signals” or “parameters” in the filename
- **Servo-U:** .sta files
- **Mindray:** .xml or .csv files
- **Syringe:** Files containing “syringe” in the filename
- **Other:** .csv or .parquet files (auto-discovers columns with datetime values)

Patient folder structure example

Figure 3: Patient folder structure example

## 4 Loading Database Options

Database options define **which data sources to enable** and **how signals should be displayed** (labels, units, colors, grouping).

### 4.1 Option 1: Default Visualization (Quick Start)

Click the green **“Default visualization (all sources)”** button. This automatically enables all 9 data sources with their default display settings. No configuration file is needed.

This is the recommended starting point for new users.

Default visualization button

Figure 4: Default visualization button

### 4.2 Option 2: Custom Configuration File

Click the blue **“Upload config file”** button to load a custom `database_options.json` file. This gives you full control over which sources are enabled and how each signal is displayed.

See Section 9 for a detailed description of the configuration file format.

Upload config file button

Figure 5: Upload config file button

## 5 Configuring Patient Options

After loading database options, the **Patient Options** form appears. It is divided into two parts: global options and per-source options.

### 5.1 Global Options

These apply to all data sources:

| Option                                    | Description                                                                                                                                            |
|-------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Path to data (folder)</b>              | Full path to the patient's root data folder                                                                                                            |
| <b>Time start filter</b>                  | Start of the time window to display (format: YYYY-MM-DD HH:MM:SS). Leave empty to use all available data.                                              |
| <b>Time end filter</b>                    | End of the time window to display. Leave empty to use all available data.                                                                              |
| <b>Re-use data if already loaded once</b> | When checked, reuses previously cached <code>.parquet</code> files from the <code>tdv_visu/</code> folder, significantly speeding up subsequent loads. |

Global patient options

Figure 6: Global patient options

### 5.2 Per-Source Options

Below the global options, each enabled data source may have additional settings displayed in individual cards arranged in a two-column grid. Common per-source options include:

- **Time shift** (seconds): Adjust the time alignment of a source relative to others. Useful when devices were not perfectly synchronized.
- **Day**: Specify the recording date for sources that require it (e.g., EIT data).

Only data sources present in the loaded database options will show their configuration cards.

Per-source options

Figure 7: Per-source options



## 6 Processing the Visualization

Once you have configured the patient options, click the large orange “**Process visualization**” button to start generating the plots.

### 6.1 What Happens During Processing

1. **Validation:** The application verifies that all mandatory fields are filled in and that the data folder exists.
2. **Data Discovery:** For each enabled data source, the application scans the patient folder for matching subfolders and files.
3. **Data Loading:** Raw data files are parsed according to each source’s format.
4. **Formatting:** Signals are filtered, resampled, and converted using your database options (labels, units, time range).
5. **Caching:** Processed data is saved as `.parquet` files in the `tdv_visu/` subfolder for faster reloading next time.
6. **Plot Generation:** Interactive Plotly figures are created and displayed in the visualization area.

A success message appears when processing completes. If no data is found for a source, it is silently skipped.

Processing visualization

Figure 8: Processing visualization

## 7 Interacting with Plots

### 7.1 Navigation Controls

Each plot provides a toolbar (top-right corner) with the following tools:

| Tool                      | Action                                           |
|---------------------------|--------------------------------------------------|
| <b>Zoom</b>               | Click and drag to zoom into a rectangular region |
| <b>Pan</b>                | Click and drag to move the view                  |
| <b>Zoom In / Zoom Out</b> | Incremental zoom buttons                         |
| <b>Autoscale</b>          | Reset the view to fit all data                   |
| <b>Reset Axes</b>         | Return to the original view                      |
| <b>Download as PNG</b>    | Save the current plot view as an image           |

### 7.2 Zooming and Panning

- **Scroll wheel:** Zoom in/out on the x-axis.
- **Click and drag:** Select a region to zoom into.
- **Double-click:** Reset the axes to show all data.

### 7.3 Dynamic Resampling (FigureResampler)

For high-frequency signals (e.g., waveforms sampled at hundreds of Hz), the application uses **Plotly-Resampler** to dynamically load detail as you zoom in. When viewing a long time range, the plot shows a downsampled overview. As you zoom into a shorter time window, the full-resolution data is loaded automatically.

This keeps the interface responsive even with millions of data points.

Interactive plot navigation

Figure 9: Interactive plot navigation

## 8 Annotations and Shape Tools

The annotation system lets you mark events, time periods, or regions of interest directly on the plots. Annotations are saved to an `annotations.json` file in the `tdv_visu/` folder and persist across sessions.

### 8.1 Drawing Annotations

Use the Plotly drawing tools in each plot's toolbar:

- **Draw Line:** Click two points to draw a vertical line marking a specific event.
- **Draw Rectangle:** Click and drag to highlight a time region or value range.

Each shape can be given a name and a color via the edit popup that appears after drawing.

### 8.2 Shape Management

After processing, the **Shape Controls** section becomes visible below the Process button:

- **Shape Dropdown:** Lists all annotations across all figures, showing the figure name and shape label.
- **Modify Button:** Opens the edit popup for the selected annotation (change name, color, or whether it spans all subplots).
- **Delete Button:** Removes the selected annotation.

### 8.3 Annotation Properties

Each annotation has the following properties:

- **Name:** A text label displayed on the shape.
- **Color:** The line or fill color.
- **Global:** When enabled, the shape spans all subplots in the figure (using paper y-coordinates).

### 8.4 Persistence

Annotations are automatically saved to `annotations.json` in the patient's `tdv_visu/` folder whenever you create, modify, or delete a shape. They are reloaded when you re-process the same patient data.

Annotation tools and shape management

Figure 10: Annotation tools and shape management

## 9 Configuration File Reference

### 9.1 patient\_options.json

This file defines patient-specific settings.

```
{
  "data_folder": "/path/to/patient/data",
  "datetime_start": "2024-10-08 10:00:00",
  "datetime_end": "2024-10-08 12:00:00",
  "quick_load": false,
  "philips_waves": {
    "time_shift": 20.0
  },
  "eit": {
    "day": "2024-10-08"
  }
}
```

| Key            | Type           | Description                                    |
|----------------|----------------|------------------------------------------------|
| data_folder    | string         | Path to the patient's root data folder         |
| datetime_start | string or null | Start of the time window (YYYY-MM-DD HH:MM:SS) |
| datetime_end   | string or null | End of the time window                         |
| quick_load     | boolean        | Reuse cached parquet files                     |
| <source_name>  | object         | Per-source options (e.g., time_shift, day)     |

### 9.2 database\_options.json

This file controls which data sources are active and how signals are displayed.

#### 9.2.1 Top-Level Structure

```
{
  "global": {
    "grouped_fields": { "Group Name": ["signal1", "signal2"] }
  },
  "philips_waves": { ... },
  "philips_numerics": { ... },
  "eit": { ... }
}
```

Each data source key is optional – only include the sources you want to enable.

#### 9.2.2 Per-Source Fields

| Key                       | Description                                       |
|---------------------------|---------------------------------------------------|
| field_display             | Array of signal names to show (others are hidden) |
| data.label_correspondence | Map raw signal names to display labels            |
| data.unit_conversion      | Multiplication factors for unit conversion        |
| data.unit_info            | Display unit strings (e.g., “mmHg”, “cmH2O”)      |

| Key                                     | Description                                                                    |
|-----------------------------------------|--------------------------------------------------------------------------------|
| <code>data.unit_range</code>            | Fixed Y-axis ranges as <code>[min, max]</code>                                 |
| <code>data.color</code>                 | Custom colors per signal                                                       |
| <code>data.priority</code>              | Plot ordering priority (lower = higher on page)                                |
| <code>data.period_resampling</code>     | Resampling period in seconds                                                   |
| <code>grouped_fields</code>             | Group signals onto the same subplot                                            |
| <code>loop</code>                       | Define PV-loop plots (e.g., <code>{"pv_loop": ["Pressure", "Volume"]}</code> ) |
| <code>numerics.period_resampling</code> | Resampling for numeric parameters                                              |
| <code>numerics.priority</code>          | Plot priority for numerics                                                     |

### 9.2.3 Global Fields

| Key                                | Description                                                     |
|------------------------------------|-----------------------------------------------------------------|
| <code>global.grouped_fields</code> | Group signals from different data sources onto the same subplot |

See `example/option_files/` in the source repository for complete example files.

## 10 Troubleshooting

### 10.1 Browser Does Not Open Automatically

If the browser does not open after launching the application, manually navigate to:

`http://127.0.0.1:8050`

Ensure no other application is using port 8050. If needed, close the terminal window which was opened in the app and restart the application.

### 10.2 No Data Found

If the visualization is empty or a data source shows no signals:

- Verify that the **data folder path** is correct and accessible.
- Check that subfolders follow the **naming conventions** (see Section 3).
- Ensure the subfolder contains files in the **expected format** for that data source.
- Check that the data source is **enabled** in your database options (or use “Default visualization” to enable all).

### 10.3 Slow Loading

Large datasets may take time to load on the first run. To speed up subsequent loads:

- Enable the “**Re-use data if already loaded once**” (`quick_load`) option. This uses the cached `.parquet` files in `tdv_visu/` instead of re-reading raw data files.

### 10.4 Time Alignment Issues

If signals from different sources appear misaligned in time:

- Use the **Time shift** option in the per-source settings to adjust alignment.
- Verify that the correct **day** or **date** is set for sources that require it (e.g., EIT).

### 10.5 Application Crashes or Errors

- Check the terminal window for error messages.
- Log files are available in the `logs/` directory (if running from source).
- Ensure the data files are not corrupted or truncated.

## 11 Appendix: Supported Data Sources

| Source             | Module Name                     | Keywords                         | File Types                  | Typical Signals                               |
|--------------------|---------------------------------|----------------------------------|-----------------------------|-----------------------------------------------|
| Philips Waves      | <code>philips_waves</code>      | <code>philips, waves</code>      | <code>.parquet</code>       | ART, PAP, CO2, respiratory pressure/volume    |
| Philips Numerics   | <code>philips_numerics</code>   | <code>philips, numerics</code>   | numerics in name            | Heart rate, SpO2, FiO2, blood pressure        |
| EIT                | <code>eit</code>                | <code>eit</code>                 | <code>.asc</code>           | Global/local impedance, impedance percentages |
| FluxMed Signals    | <code>fluxmed_signals</code>    | <code>fluxmed, signals</code>    | signals in name             | Respiratory waveforms                         |
| FluxMed Parameters | <code>fluxmed_parameters</code> | <code>fluxmed, parameters</code> | parameters in name          | Respiratory parameters                        |
| Servo-U            | <code>servo_u</code>            | <code>servo</code>               | <code>.sta</code>           | Ventilator waveforms and settings             |
| Mindray            | <code>mindray</code>            | <code>mindray</code>             | <code>.xml, .csv</code>     | Monitor waveforms and parameters              |
| Syringe            | <code>syringe</code>            | <code>syringe</code>             | syringe in name             | Infusion rates and volumes                    |
| Other              | <code>other</code>              | <code>other</code>               | <code>.csv, .parquet</code> | Any time-series with datetime columns         |